

Data Sheet:

WanderGuard[®] BLUE Tag

Signals the Controller when an at-risk resident approaches and can lock a monitored door

Description

The WanderGuard BLUE Tag signals the WanderGuard BLUE EX5700 Controller that it is in field. The Controller issues an alarm and can lock a monitored door so that an at-risk resident wearing the Tag cannot leave through the door.

The WanderGuard BLUE Tag is a component of STANLEY Healthcare's Wanderguard BLUE Wander Management Solution for monitoring facility doors and preventing at-risk resident wandering. The resident typically wears the WanderGuard BLUE Tag on the wrist like a bracelet.

How It Works

Tags are activated by the WanderGuard BLUE Detector tool. When an at-risk resident wearing the Tag approaches the door and is in the LF range of the WanderGuard BLUE Controller, the Tag responds and transmits BLE signals so that the Controller can immediately lock the door.

Using the WanderGuard BLUE Detector, information about the Tag including Activation date, Battery level, Status and Tag ID can be displayed in the WanderGuard BLUE Manager mobile application.

The Tag ID is printed on the side of the Tag.









Product Highlights

- Each Tag is unique Unique electronic serial number to identify each Tag.
- Wireless Activated wirelessly by the WanderGuard BLUE Detector.
- Tag data visible in the mobile application – Tag data can be viewed in WanderGuard BLUE Manager.
- Ninety-day or three-year Tags.
- Multiple Tag activation Multiple Tags can be simultaneously activated.
- Attractive, high-impact, tamperresistant case.
- Standard and rugged wristbands Standard wrist bands are provided. Resilient and comfortable, cutproof Securaband[™] wrist bands are optional.
- Easily sterilized Smooth surface permits easy sterilization and prevents harboring of bacteria.
- Waterproof IP67-compliant.
- Small and unobtrusive One of the smallest and lightest tags on the market. Worn like a bracelet.



Product Specifications	
Part Number	SKU for Ninety-day Tag: WGB-TAG-1000-90D
	SKU for Three-year Tag: WGB-TAG-1000-3Y
Model	WGB-TAG-1000-90D
	WGB-TAG-1000-3Y
LF Frequency	125 KHz
BLE Frequency Range	2400-2483.5 MHz
Battery Life Options	Ninety days
	Three years
Temperature	32° to 122°F (0 to 50°C)
Humidity	0-95% RH @ 70°F (21°C), non-condensing
Compliance	IP67
Dimensions	0-95% RH @ 70°F (21°C), non-condensing
Weight	Approximately 0.41 oz. (11.6 g)

By Your Side™ Lifetime Customer Care

At STANLEY Healthcare, we measure our success by yours. The By Your Side Customer Care commitment combines a range of services to ensure that every customer achieves success and realizes the full value of their investment.



STANLEY Healthcare 130 Turner Street, Waltham, MA 02453 Phone: +1-888-622-6992 Email: stanleyhealthcare@sbdinc.com

© 2017 STANLEY Healthcare. All rights reserved. 0971-291-000 Rev A. 2017/07/20. KB Article: 10158.